

DyLight® 650 Anti-T7 tag® antibody ab119353

製品の概要

製品名	DyLight® 650 Anti-T7 tag® antibody
製品の詳細	DyLight® 650 Rabbit polyclonal to T7 tag®
由来種	Rabbit
標識	DyLight® 650. Ex: 654nm, Em: 673nm
アプリケーション	適用あり: ICC/IF
種交差性	交差種: Species independent
免疫原	Synthetic peptide: MASMTGGQMG corresponding to T7 tag®.

 [Run BLAST with](#)

 [Run BLAST with](#)

特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. The product can be stored for up to 12 months.
バッファー	pH: 8.20 Constituent: 99% Tris buffered saline
精製度	Protein A purified
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab119353の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/500 - 1/5000.

ターゲット情報

関連性

The T7 tag is an 11 amino acid peptide encoded in the leader sequence of T7 bacteriophage gene10. This gene encodes a T7 major capsid protein whose function is not clear. The T7 tag serves as a tag in many expression vectors including the pET system that is based on the very efficient T7 RNA polymerase expression system. Monoclonal antibodies specific for T7 tag are an important tool for studying expression of recombinant T7-tagged proteins.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors