

# Anti-EBV Nuclear Antigen/EBNA1 antibody ab20870

### 製品の概要

製品名	Anti-EBV Nuclear Antigen/EBNA1 antibody
製品の詳細	Goat polyclonal to EBV Nuclear Antigen/EBNA1
由来種	Goat
特異性	Epstein-Barr nuclear antigen. Non reactive with human lymphoid cells. Functions in ELISA against various recombinant EBNA-1 preps.
アプリケーション	<b>適用あり:</b> ELISA, ICC/IF
種交差性	<b>交差種:</b> Epstein-Barr virus
免疫原	Recombinant full length protein corresponding to EBV Nuclear Antigen/EBNA1.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
精製度	Ion Exchange Chromatography
特記事項(精製)	>95% pure. Sodium sulfate precipitation & ion-exchange chromatography.
ポリモノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab20870の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.

## ターゲット情報

関連性	Epstein-Barr virus (EBV) nuclear antigen 1 (EBNA1) is the one EBV antigen that is expressed in all EBV associated malignancies. It has long been thought to go undetected by the cell mediated immune system. However, recent studies show that EBNA1 can be presented to both CD4+ and CD8+ T cells, making it a potential new target for immunotherapy of EBV related cancers.
細胞内局在	Nuclear. Free in the nucleoplasm, somewhat associated with the chromatin and hardly, if at all associated with the nuclear matrix.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors