

Anti-Eph receptor B4/HTK antibody [RM0056-6A14] ab73259

★★★★★ [1 Abreviews](#) [4 References](#) [画像数 2](#)

製品の概要

製品名	Anti-Eph receptor B4/HTK antibody [RM0056-6A14]
製品の詳細	Rat monoclonal [RM0056-6A14] to Eph receptor B4/HTK
由来種	Rat
特異性	ab73259 is specific for mouse Eph receptor B4/HTK. It does not cross-react with recombinant mouse EphA2, EphA3, EphA4, EphA6, EphA7 and EphB6
アプリケーション	適用あり: IHC-P
種交差性	交差種: Mouse
免疫原	Recombinant fragment corresponding to Mouse Eph receptor B4/HTK (extracellular). Database link: P54761
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituent: PBS
精製度	Protein A purified
特記事項 (精製)	ab73259 was produced from a hybridoma (mouse myeloma fused with spleen cells from a rat immunized with recombinant mouse Eph receptor B4/HTK extracellular domain).
ポリモノ	モノクローナル
クローン名	RM0056-6A14
アイソタイプ	IgG1

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab73259の使用に適用されず**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (1)	1/100 - 1/200. Perform enzymatic antigen retrieval with 20 ug/ml proteinase K before commencing with IHC protocol.

ターゲット情報


機能	Receptor for members of the ephrin-B family. Binds to ephrin-B2. May have a role in events mediating differentiation and development.
組織特異性	Abundantly expressed in placenta and in a range of primary tissues and malignant cell lines. Expressed in fetal, but not adult, brain, and in primitive and myeloid, but not lymphoid, hematopoietic cells.
配列類似性	Belongs to the protein kinase superfamily. Tyr protein kinase family. Ephrin receptor subfamily. Contains 2 fibronectin type-III domains. Contains 1 protein kinase domain. Contains 1 SAM (sterile alpha motif) domain.
翻訳後修飾	Autophosphorylated.
細胞内局在	Membrane.

画像



ab73259 at 1/200 dilution, staining Eph receptor B4/HTK in LPS treated mouse kidney tissue section by Immunohistochemistry (Formalin/PFA fixed paraffin-embedded sections).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Eph receptor B4/HTK antibody [RM0056-6A14] (ab73259)



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Eph receptor B4/HTK antibody [RM0056-6A14] (ab73259)

This image is courtesy of an anonymous Abreview

ab73259 staining the Eph receptor B4/HTK in Mouse lymph node tissue sections by Immunohistochemistry (IHC-P - paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with paraformaldehyde and blocked with 10% serum for 1 hour at 20°C; antigen retrieval was not performed. Samples were incubated with primary antibody (1/100) for 12 hours at 4°C. A Biotin-conjugated Rabbit anti-rat polyclonal (1/200) was used as the secondary antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors