

# Anti-ESAT6 antibody ab45073

★☆☆☆☆ **3 Abreviews** **4 References**

### 製品の概要

製品名	Anti-ESAT6 antibody
製品の詳細	Rabbit polyclonal to ESAT6
由来種	Rabbit
特異性	The rabbit polyclonal antibody to EsaT-6 recognizes the EsaT-6 (Rv3875) protein of Mycobacterium tuberculosis and Mycobacterium bovis.
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Mycobacterium bovis, Mycobacterium tuberculosis
免疫原	Recombinant full length protein corresponding to ESAT6.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituent: PBS</p>
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab45073の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★☆☆☆☆ (2)	Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 11 kDa. Use under reducing conditions.

## ターゲット情報

関連性	The ESAT6 antigen from Mycobacterium tuberculosis is a dominant target for cell mediated immunity in the early phase of tuberculosis (TB) in TB patients as well as in various animal models. ESAT6 is not found in M. bovis BCG. Its function is not known but it elicits high level of IFN-gamma from memory effector cells during the first phase of a protective immune response.
細胞内局在	Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors