

# Anti-Ethinylestradiol antibody ab59670

## 1 References

### 製品の概要

製品名	Anti-Ethinylestradiol antibody
製品の詳細	Sheep polyclonal to Ethinylestradiol
由来種	Sheep
アプリケーション	<b>適用あり:</b> ELISA, EIA
種交差性	<b>交差種:</b> Species independent
免疫原	Ethinylestradiol conjugated to Bovine thyroglobulin
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: 0.0536% PBS, 0.87% Sodium chloride</p>
精製度	IgG fraction
特記事項 (精製)	Ig fraction prepared by caprylic acid and ammonium sulphate precipitation
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab59670の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use a concentration of 1.25 µg/ml.
EIA		Use at an assay dependent concentration.

## ターゲット情報

関連性	Ethinylestradiol, (8R,9S,13S,14S,17R)-17-ethynyl-13-methyl-7,8,9,11,12,14,15,16-octahydro-6H-cyclopenta[a]phenanthrene-3,17-diol, is a semisynthetic alkylated Estradiol with a 17-alpha-ethynyl substitution. It has high estrogenic potency when administered orally, and is often used as the estrogenic component in oral contraceptives. It has a molecular weight 296.40336 g/mol and molecular formula C <sub>20</sub> H <sub>24</sub> O <sub>2</sub> .
-----	--

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors