

Anti-Exoenzyme S antibody ab20031

製品の概要

製品名	Anti-Exoenzyme S antibody
製品の詳細	Chicken polyclonal to Exoenzyme S
由来種	Chicken
アプリケーション	適用あり: WB
種交差性	交差種: Pseudomonas aeruginosa
免疫原	Fusion protein, corresponding to amino acids 366-453 of Exoenzyme S strain PAO1
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Sodium azide
精製度	Protein L purified
特記事項(精製)	Total IgY.
ポリ/モノ	ポリクローナル
アイソタイプ	IgY

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab20031の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/5000.

ターゲット情報

関連性	Pseudomonas aeruginosa produces two extracellular protein toxins, Exoenzyme S and Exotoxin A. Exoenzyme S is probably an exotoxin. It has the characteristic subunit structure of the A-component of a bacterial toxin, and it has ADP-ribosylating activity (for a variety of eukaryotic proteins) characteristic of exotoxins. It has been suggested that exoenzyme S may act to impair the function of phagocytic cells in the bloodstream and internal organs to prepare for invasion by Ps. aeruginosa. Exoenzyme S is a potent stimulus for the proliferation of T cells in greater than 95% of adults, which could contribute to the pathogenesis of Cystic Fibrosis.
細胞内局在	Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors