

Anti-F4/80 antibody [Cl:A3-1] - Macrophage Marker ab6640

★★★★★ [69 Abreviews](#) [961 References](#) [画像数 5](#)

製品の概要

製品名	Anti-F4/80 antibody [Cl:A3-1] - Macrophage Marker
製品の詳細	Rat monoclonal [Cl:A3-1] to F4/80 - Macrophage Marker
由来種	Rat
特異性	Anti-F4/80 antibody (ab6640) recognises the mouse F4/80 antigen, a 160kD glycoprotein expressed by murine macrophages.
アプリケーション	適用あり: Flow Cyt, ICC/IF
種交差性	交差種: Mouse
免疫原	Tissue, cells or virus corresponding to F4/80. Thioglycollate stimulated peritoneal macrophages from C57/BL mice
ポジティブ・コントロール	ICC/IF: Raw 264.7 cells, mouse bone marrow cells. Flow Cyt: J774 cells and peritoneal macrophages.
特記事項	<p>Please note we cannot guarantee IHC-P, however should you use this application for this product then please use a specific protocol which can be found in the protocols section of our datasheet here.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Protein G purified
特記事項 (精製)	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant.

ポリ/モノ	モノクローナル
クローン名	Cl:A3-1
ミエローマ	NS1
アイソタイプ	IgG2b

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab6640の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		1/10 - 1/50. ab18536 - Rat monoclonal IgG2b, is suitable for use as an isotype control with this antibody.
ICC/IF	★★★★★ (11)	Use a concentration of 5 - 10 µg/ml. We recommend ab150165 Goat Anti-Rat IgG H&L (Alexa Fluor® 488) secondary antibody.

ターゲット情報


機能	Orphan receptor involved in cell adhesion and probably in cell-cell interactions specifically involving cells of the immune system. May play a role in regulatory T-cells (Treg) development.
組織特異性	Expression is restricted to eosinophils.
配列類似性	Belongs to the G-protein coupled receptor 2 family. Adhesion G-protein coupled receptor (ADGR) subfamily. Contains 6 EGF-like domains. Contains 1 GPS domain.
細胞内局在	Cell membrane.

画像



Immunocytochemistry/ Immunofluorescence - Anti-F4/80 antibody [Cl:A3-1] - Macrophage Marker (ab6640)

ab6640 stained Raw 264.7 cells. The cells were 100% methanol fixed for 5 minutes at room temperature and then incubated in 1% BSA / 10% normal Goat serum / 0.3M glycine in 0.1% PBS-Tween for 1 hour at room temperature to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody (ab6640 at 5 µg/ml) overnight at +4°C. The secondary antibody (pseudo-colored green) was **ab150165 Goat Anti-Rat IgG H&L (Alexa Fluor® 488)** preadsorbed used at a 1/1000 dilution for 1 hour at room temperature. Alexa Fluor® 594 WGA was used to label plasma membranes (pseudo-colored red) at a 1/200 dilution for 1 hour at room temperature. DAPI was used to stain the cell nuclei (pseudo-colored blue) at a concentration of 1.43 µM for 1 hour at room temperature.



Flow Cytometry - Anti-F4/80 antibody [Cl:A3-1] - Macrophage Marker (ab6640)


Flow cytometry analysis of J774 cells labelling F4/80 (red) with ab6640 at a dilution of 1/10 followed by streptavidin FITC secondary at 1/100 dilution. The blue line shows J774 cells stained with Rat anti mouse isotype control.



Immunocytochemistry/ Immunofluorescence - Anti-F4/80 antibody [Cl:A3-1] - Macrophage Marker (ab6640)


This image is courtesy of an anonymous Abreview

ab6640 staining F4/80 in mouse bone marrow cells by immunocytochemistry/ immunofluorescence. Cells were formaldehyde fixed and permeabilized in 0.2% Triton X-100 prior to blocking in 2% BSA for 30 minutes at 20°C. The primary antibody was diluted 1/200 and incubated with the sample for 9 hours at 4°C. Alexa fluor® 488 goat polyclonal to rat Ig, diluted 1/200, was used as the secondary antibody.



Flow Cytometry - Anti-F4/80 antibody [Cl:A3-1] -
Macrophage Marker (ab6640)

Flow cytometry analysis of peritoneal macrophages labelling F4/80 (red) with ab6640 at a dilution of 1/10. The blue line shows J774 cells stained with Rat anti mouse isotype control.



Flow Cytometry - Anti-F4/80 antibody [Cl:A3-1] -
Macrophage Marker (ab6640)

Flow cytometry analysis of J774 cells labelling F4/80 (red) with ab6640 at a dilution of 1/50 followed by goat anti rat IgG FITC secondary antibody at 1/100 dilution. The blue line shows J774 cells stained with Rat anti mouse isotype control.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors