

Anti-FCN1/M-Ficolin antibody ab70814

画像 1

医薬用外劇物

製品の概要

製品名	Anti-FCN1/M-Ficolin antibody
製品の詳細	Rabbit polyclonal to FCN1/M-Ficolin
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Rat
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	3T3 cell lysate and rat kidney tissue lysate
特記事項	<p>This product is manufactured by BioVision, an Abcam company and was previously called 3852 Ficolin-1 Antibody. 3852-100 is the same size as the 100 µg size of ab70814.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.2 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.5% BSA, 30% Glycerol (glycerin, glycerine)
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab70814の使用に適用されず
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.5 - 4 µg/ml. Detects a band of approximately 33 , 37 kDa (predicted molecular weight: 35 kDa).

ターゲット情報

機能	Involved in serum exerting lectin activity. Binds GlcNAc.
組織特異性	Peripheral blood leukocytes. Also detected in spleen, lung, and thymus, may be due to the presence of tissue macrophages or trapped blood in these tissues. Not detected on lymphocytes and granulocytes.
配列類似性	Belongs to the ficolin lectin family. Contains 1 collagen-like domain. Contains 1 fibrinogen C-terminal domain.
細胞内局在	Secreted. Found on the monocyte surface.

画像



All lanes : Anti-FCN1/M-Ficolin antibody (ab70814) at 0.5 µg/ml

Lane 1 : 3T3 cell lysate

Lane 2 : Rat kidney tissue lysate

Predicted band size: 35 kDa

Observed band size: 33,37 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors