

Anti-Feline Leukemia Virus (g27) antibody [YCH] ab116669

2 References

製品の概要

製品名	Anti-Feline Leukemia Virus (g27) antibody [YCH]
製品の詳細	Mouse monoclonal [YCH] to Feline Leukemia Virus (g27)
由来種	Mouse
特異性	ab116669 targets Feline Leukemia Virus (g27) in the tested applications and shows reactivity with Virus samples.
アプリケーション	適用あり: IP, ELISA, RIA
種交差性	交差種: Feline leukemia virus
免疫原	Tissue, cells or virus corresponding to Feline Leukemia Virus (g27).
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	pH: 7.4 Preservative: 0.05% Sodium azide Constituents: 0.019% Potassium chloride, 0.79% Sodium chloride, 0.135% Potassium phosphate, 0.1% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	YCH
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab116669の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IP		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.
RIA		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors