

### Anti-FGL2/Prothrombinase antibody ab198029

★☆☆☆☆ [1 Abreviews](#) [4 References](#) [画像数 2](#)

#### 製品の概要

製品名	Anti-FGL2/Prothrombinase antibody
製品の詳細	Rabbit polyclonal to FGL2/Prothrombinase
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IHC-P, WB
種交差性	<b>交差種:</b> Mouse, Human
免疫原	Recombinant fragment (GST-tag) corresponding to Human FGL2/Prothrombinase (internal sequence). BC033820 Database link: <a href="#">Q14314</a>
ポジティブ・コントロール	Mouse thymus tissue lysate; Human liver cancer tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 7.4 Preservative: 0.05% Sodium azide Constituents: 50% Glycerol (glycerin, glycerine), 49% PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

## The Abpromise guarantee


**Abpromise保証は、**次のテスト済みアプリケーションにおけるab198029の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/20 - 1/100.
WB		1/450 - 1/2000. Predicted molecular weight: 50 kDa.

## ターゲット情報

機能	May play a role in physiologic lymphocyte functions at mucosal sites.
組織特異性	Constitutively expressed in cytotoxic T-cells.
配列類似性	Contains 1 fibrinogen C-terminal domain.
細胞内局在	Secreted.

## 画像




Anti-FGL2/Prothrombinase antibody (ab198029) at 1/450 dilution +  
Mouse thymus tissue lysate at 40 µg

**Predicted band size:** 50 kDa

**Exposure time:** 40 seconds

8% SDS-PAGE



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue labeling FGL2/Prothrombinase with ab198029 at 1/20 dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors