

FITC Anti-CD137 antibody [4B4 (4B4-1)] ab93520

1 References [画像数 1](#)

製品の概要

製品名	FITC Anti-CD137 antibody [4B4 (4B4-1)]
製品の詳細	FITC Mouse monoclonal [4B4 (4B4-1)] to CD137
由来種	Mouse
標識	FITC. Ex: 493nm, Em: 528nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	Tissue, cells or virus corresponding to Human CD137. Isolated Human eosinophils.
ポジティブ・コントロール	ConA-stimulated normal Human peripheral blood cells
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.20
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	4B4 (4B4-1)
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab93520の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 5µl for 10 ⁵⁻⁸ cells. ab106163 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody. Use 5µl for 10 ⁵⁻⁸ cells in a final volume of 100µl.

ターゲット情報

機能	Receptor for TNFSF14/4-1BBL. Possibly active during T cell activation.
組織特異性	Expressed on the surface of activated T-cells.
配列類似性	Contains 4 TNFR-Cys repeats.
細胞内局在	Membrane.

画像



Staining of ConA-stimulated human peripheral blood cells with Anti-Human CD69 APC and Mouse IgG1 kappa Isotype Control FITC (left) or ab93520 (right). Total viable cells were used for analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors