

### FITC Anti-CD28 antibody [PV-1] ab24844

画像数 1

#### 製品の概要

製品名	FITC Anti-CD28 antibody [PV-1]
製品の詳細	FITC Armenian hamster monoclonal [PV-1] to CD28
由来種	Armenian hamster
標識	FITC. Ex: 493nm, Em: 528nm
アプリケーション	<b>適用あり:</b> Flow Cyt
種交差性	<b>交差種:</b> Mouse, Human
免疫原	Tissue, cells or virus corresponding to Mouse CD28. C57BL/6N mouse T cell lymphoma EL4 cell line
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	<p>pH: 7.4</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituent: PBS</p>
精製度	Affinity purified
ポリ/モノ	モノクローナル
クローン名	PV-1
アイソタイプ	IgG

#### アプリケーション

**The Abpromise guarantee**      Abpromise保証は、次のテスト済みアプリケーションにおけるab24844の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 3µg for 10 <sup>6</sup> cells. <b>ab18473</b> - Armenian Hamster monoclonal IgG, is suitable for use as an isotype control with this antibody.

## ターゲット情報

機能	Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.
組織特異性	Expressed in T-cells and plasma cells, but not in less mature B-cells.
配列類似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
細胞内局在	Membrane.

## 画像



ab24844 at a concentration of 3 µg/10<sup>6</sup> cells staining human A 20 (B cell lymphoma cells) and EL-4 cells by flow cytometry. A 50:50 mixture of the A20 B-cell lymphoma and the EL-4 T-cell lymphoma was stained with ab24844. Large cells were gated and analyzed by flow cytometry.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors