

# FITC Anti-CD5 antibody [M28623] ab52310

## 1 References

### 製品の概要

製品名	FITC Anti-CD5 antibody [M28623]
製品の詳細	FITC Mouse monoclonal [M28623] to CD5
由来種	Mouse
標識	FITC. Ex: 493nm, Em: 528nm
アプリケーション	<b>適用あり:</b> ICC/IF, Flow Cyt
種交差性	<b>交差種:</b> Human
免疫原	Tissue, cells or virus corresponding to Human CD5. Database link: <a href="#">P06127</a>
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.08% Sodium azide Constituent: PBS
	0.2% protein carrier
ポリ/モノ	モノクローナル
クローン名	M28623
アイソタイプ	IgG2a
軽鎖の種類	kappa

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab52310の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent concentration.
Flow Cyt		Use 10µl for 10 <sup>6</sup> cells. <b>ab18455</b> - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody. Use 10µl for 10 <sup>6</sup> cells or 100µl whole blood.

## ターゲット情報

機能	May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.
配列類似性	Contains 3 SRCR domains.
細胞内局在	Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors