

FITC Anti-Cytomegalovirus Glycoprotein B antibody [HCMV37] ab74085

1 References

製品の概要

製品名	FITC Anti-Cytomegalovirus Glycoprotein B antibody [HCMV37]
製品の詳細	FITC Mouse monoclonal [HCMV37] to Cytomegalovirus Glycoprotein B
由来種	Mouse
標識	FITC. Ex: 493nm, Em: 528nm
特異性	ab74085 reacts with MRC-5 infected cells but not uninfected cells. This antibody has been found to be neutralizing.
アプリケーション	適用あり: ICC/IF
種交差性	交差種: Human cytomegalovirus
免疫原	Tissue, cells or virus corresponding to Cytomegalovirus Glycoprotein B. Human CMV extracts.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: 1% BSA, PBS</p>
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	HCMV37
アイソタイプ	IgG2a

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab74085の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent dilution.

ターゲット情報

関連性 Human cytomegalovirus (HCMV) glycoprotein B (gB) is an abundant virion envelope protein that has been shown to be essential for the infectivity of HCMV. HCMV gB is also one of the most immunogenic virus-encoded proteins, and a significant fraction of virus neutralizing antibodies are directed at gB.

細胞内局在 Virus envelope.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors