


# FITC Anti-DMPO Nitron Adduct antibody [N1664A] ab110429

## 1 References

### 製品の概要

製品名	FITC Anti-DMPO Nitron Adduct antibody [N1664A]
製品の詳細	FITC Mouse monoclonal [N1664A] to DMPO Nitron Adduct
由来種	Mouse
標識	FITC. Ex: 493nm, Em: 528nm
特異性	ab110429 recognizes DMPO, DMPO-octanoic acid, DMPO-protein adducts and DMPO-DNA adducts. Does not cross-react with non-adducted proteins or DNA.
アプリケーション	適用あり: ICC
種交差性	<p>交差が予測される動物種: Mouse </p>
免疫原	Chemical/ Small Molecule corresponding to DMPO Nitron Adduct conjugated to ovalbumin.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.20 Constituents: 0.27% Potassium phosphate, 0.87% Sodium chloride, 0.1% BSA
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	N1664A
アイソタイプ	IgG1

## アプリケーション

**The Abpromise guarantee**      Abpromise保証は、次のテスト済みアプリケーションにおけるab110429の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC		1/100.

## ターゲット情報

**関連性**      For measurement of protein radicals; some proteins with endogenous peroxidase activity are especially susceptible to forming DMPO adducts.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors