

### FITC Anti-IL-2 antibody [JES6-5H4] ab87029

#### 製品の概要

---

製品名	FITC Anti-IL-2 antibody [JES6-5H4]
製品の詳細	FITC Rat monoclonal [JES6-5H4] to IL-2
由来種	Rat
標識	FITC. Ex: 493nm, Em: 528nm
アプリケーション	<b>適用あり:</b> ELISA, Flow Cyt
種交差性	<b>交差種:</b> Mouse
免疫原	Recombinant full length protein corresponding to Mouse IL-2.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: PBS
	Also contains a stabilizing agent.
精製度	IgG fraction
ポリ/モノ	モノクローナル
クローン名	JES6-5H4
アイソタイプ	IgG2b
軽鎖の種類	kappa

**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab87029の使用に適用されず**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration.

## ターゲット情報

<b>機能</b>	Produced by T-cells in response to antigenic or mitogenic stimulation, this protein is required for T-cell proliferation and other activities crucial to regulation of the immune response. Can stimulate B-cells, monocytes, lymphokine-activated killer cells, natural killer cells, and glioma cells.
<b>関連疾患</b>	Note=A chromosomal aberration involving IL2 is found in a form of T-cell acute lymphoblastic leukemia (T-ALL). Translocation t(4;16)(q26;p13) with involves TNFRSF17.
<b>配列類似性</b>	Belongs to the IL-2 family.
<b>細胞内局在</b>	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors