

### FITC Anti-IL-2 Receptor alpha antibody [3C7] ab25001

#### 製品の概要

製品名	FITC Anti-IL-2 Receptor alpha antibody [3C7]
製品の詳細	FITC Rat monoclonal [3C7] to IL-2 Receptor alpha
由来種	Rat
標識	FITC. Ex: 493nm, Em: 528nm
特異性	This antibody is specific to the low affinity alpha chain IL-2 receptor.
免疫原	Tissue, cells or virus corresponding to IL-2 Receptor alpha. IL-2 dependent BALB/c helper T cell clone HT2
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	<p>pH: 7.4</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituent: PBS</p>
精製度	Affinity purified
ポリ/モノ	モノクローナル
クローン名	3C7
アイソタイプ	IgG2b
軽鎖の種類	kappa

#### ターゲット情報

機能	Receptor for interleukin-2.
----	-----------------------------

関連疾患	Genetic variations in IL2RA are associated with susceptibility to diabetes mellitus insulin-dependent type 10 (IDDM10) [MIM:601942]. A multifactorial disorder of glucose homeostasis that is characterized by susceptibility to ketoacidosis in the absence of insulin therapy. Clinical features are polydipsia, polyphagia and polyuria which result from hyperglycemia-induced osmotic diuresis and secondary thirst. These derangements result in long-term complications that affect the eyes, kidneys, nerves, and blood vessels.
配列類似性	Contains 2 Sushi (CCP/SCR) domains.
細胞内局在	Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors