

FITC Anti-Klebsiella spp antibody ab69467

製品の概要

製品名	FITC Anti-Klebsiella spp antibody
製品の詳細	FITC Rabbit polyclonal to Klebsiella spp
由来種	Rabbit
標識	FITC. Ex: 493nm, Em: 528nm
特異性	ab69467 reacts with Klebsiella sp. in bacterial and infected tissue samples.
アプリケーション	適用あり: IHC-Fr, ICC/IF
種交差性	交差種: Klebsiella pneumoniae
免疫原	Tissue, cells or virus corresponding to Klebsiella spp. Whole cells Klebsiella pneumoniae
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.1% Sodium azide Constituents: 1% BSA, PBS
精製度	Affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab69467の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-Fr		1/10 - 1/50. Fix with Acetone.
ICC/IF		1/10.

ターゲット情報

関連性	<p>Klebsiella is a genus of non-motile, Gram-negative, Oxidase-negative bacteria with a prominent polysaccharide-based capsule. Frequent human pathogens, Klebsiella organisms can lead to a wide range of disease states, notably pneumonia, urinary tract infections, septicemia, Ankylosing spondylitis, and soft tissue infections. Klebsiella pneumoniae is clinically the most important member of the Klebsiella genus of Enterobacteriaceae. New antibiotic resistant strains of K. pneumoniae are appearing, and it is increasingly found as a nosocomial infection.</p>
-----	---

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors