

FITC Mouse IgG2a, kappa monoclonal [MOPC-173] - isotype control ab18455

1 References

製品の概要

| | |
|-----------|---|
| 製品名 | FITC Mouse IgG2a, kappa monoclonal [MOPC-173] - isotype control |
| 標識 | FITC. Ex: 493nm, Em: 528nm |
| 特記事項 (標識) | conjugated with FITC under optimal conditions. |
| アプリケーション | 適用あり: Flow Cyt |
| 免疫原 | This information is considered to be commercially sensitive. |
| 特記事項 | <p>The ab18455, MOPC-173 immunoglobulin has unknown specificity and was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed rat and human tissues.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|------------------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | <p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituent: PBS</p> |
| 精製度 | Protein A purified |
| 特記事項 (精製) | The ab18455, MOPC-173 antibody was purified by affinity chromatography The solution is free of unconjugated FITC. |
| アイソタイプ・コントロール 備考 | The ab18455, MOPC-173 immunoglobulin has unknown specificity and was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed rat and human tissues. |
| ポリ/モノ | モノクローナル |

| | |
|--------|----------|
| クローン名 | MOPC-173 |
| アイソタイプ | IgG2a |
| 軽鎖の種類 | kappa |

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab18455の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| Flow Cyt | | Use at an assay dependent concentration. Suitable for surface or intracellular staining. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors