




### Anti-FOXA1 antibody ab5089

★★★★★ [2 Abreviews](#) [80 References](#) [画像数 4](#)

#### 製品の概要

製品名	Anti-FOXA1 antibody
製品の詳細	Goat polyclonal to FOXA1
由来種	Goat
アプリケーション	<b>適用あり:</b> WB, ICC
種交差性	<b>交差種:</b> Mouse, Human <b>交差が予測される動物種:</b> Rat, Cow, Pig, Xenopus laevis, Zebrafish 
免疫原	Synthetic peptide corresponding to Human FOXA1 aa 462-472 (C terminal). Sequence: C-GVYSRPVLNTS  Database link: <a href="#">P55317</a> (Peptide available as <a href="#">ab23074</a> )   <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a>
ポジティブ・コントロール	WB: HepG2 and MCF7 nuclear cell lysates; Mouse liver lysate. ICC: MCF7 and U2OS cells.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.3 Preservative: 0.02% Sodium azide Constituents: Tris buffered saline, 0.5% BSA
精製度	Ammonium Sulphate Precipitation

特記事項 (精製)	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee      **Abpromise保証は、次のテスト済みアプリケーションにおけるab5089の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	Use a concentration of 0.1 - 0.3 µg/ml. Detects a band of approximately 50 kDa (predicted molecular weight: 49.1 kDa). A 1 hour primary incubation at room temperature is recommended for this product. An additional band of unknown identity was also consistently observed at 37kDa. This band was successfully blocked by incubation with the immunizing peptide.
ICC		Use a concentration of 10 µg/ml.

ターゲット情報


機能	Transcription factor that is involved in embryonic development, establishment of tissue-specific gene expression and regulation of gene expression in differentiated tissues. Is thought to act as a 'pioneer' factor opening the compacted chromatin for other proteins through interactions with nucleosomal core histones and thereby replacing linker histones at target enhancer and/or promoter sites. Binds DNA with the consensus sequence 5'-[AC]A[AT]T[AG]TT[GT][AG][CT]T[CT]-3' (By similarity). Proposed to play a role in translating the epigenetic signatures into cell type-specific enhancer-driven transcriptional programs. Its differential recruitment to chromatin is dependent on distribution of histone H3 methylated at 'Lys-5' (H3K4me2) in estrogen-regulated genes. Involved in the development of multiple endoderm-derived organ systems such as liver, pancreas, lung and prostate; FOXA1 and FOXA2 seem to have at least in part redundant roles (By similarity). Modulates the transcriptional activity of nuclear hormone receptors. Is involved in ESR1-mediated transcription; required for ESR1 binding to the NKX2-1 promoter in breast cancer cells; binds to the RPRM promoter and is required for the estrogen-induced repression of RPRM. Involved in regulation of apoptosis by inhibiting the expression of BCL2. Involved in cell cycle regulation by activating expression of CDKN1B, alone or in conjunction with BRCA1. Originally described as a transcription activator for a number of liver genes such as AFP, albumin, tyrosine aminotransferase, PEPCK, etc. Interacts with the cis-acting regulatory regions of these genes. Involved in glucose homeostasis.
組織特異性	Highly expressed in prostate and ESR1-positive breast tumors. Overexpressed in esophageal and lung adenocarcinomas.
配列類似性	Contains 1 fork-head DNA-binding domain.
細胞内局在	Nucleus.

画像



Immunofluorescent analysis of paraformaldehyde-fixed, 0.15% TritonX-100 permeabilized MCF7 (Human breast adenocarcinoma cell line) cells labeling FOXA1 with ab5089 at 10 ug/ml (1h). Image showing nuclear staining in MCF7 cells. An Alexa Fluor 488<sup>®</sup> secondary antibody was used at 2 ug/ml. Actin filaments were stained with phalloidin (red) and the nuclear stain is DAPI (blue).

Secondary antibody only control: Unimmunized goat IgG (10 ug/ml) followed by an Alexa Fluor 488<sup>®</sup> secondary antibody at 2ug/ml.




**Lane 1 :** Anti-FOXA1 antibody (ab5089) at 0.1 µg/ml  
**Lane 2 :** Anti-FOXA1 antibody (ab5089) at 0.3 µg/ml

**Lane 1 :** HepG2 (Human liver hepatocellular carcinoma cell line) nuclear cell lysate  
**Lane 2 :** MCF7 (Human breast adenocarcinoma cell line) nuclear cell lysate

Lysates/proteins at 35 µg per lane.

**Predicted band size:** 49.1 kDa


Western blot analysis labeling FOXA1 with ab5089.



Anti-FOXA1 antibody (ab5089) at 0.1 µg/ml + Mouse liver lysate at 35 µg

**Predicted band size:** 49.1 kDa

Western blot labeling FOXA1 with ab5089 at 0.1 ug/ml.



Immunofluorescent analysis of paraformaldehyde-fixed, 0.15% TritonX-100 permeabilized U-2 OS (Human bone osteosarcoma epithelial cell line) cells labeling FOXA1 with ab5089 at 10 ug/ml (1h). Image showing nuclear staining in U-2 OS cells. An Alexa Fluor 488® secondary antibody was used at 2 ug/ml. Actin filaments were stained with phalloidin (red) and the nuclear stain is DAPI (blue).

Secondary antibody only control: Unimmunized goat IgG (10 ug/ml) followed by an Alexa Fluor 488® secondary antibody at 2ug/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors