

Anti-FOXN1 antibody ab96374

画像 1

医薬用外劇物

製品の概要

製品名	Anti-FOXN1 antibody
製品の詳細	Rabbit polyclonal to FOXN1
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to a region within amino acids 585 - 648 of Human FOXN1 (NP_003584).
ポジティブ・コントロール	HeLa, MOLT4 cell lysates and Raji whole cell lysates
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 89.99% PBS, 10% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab96374の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 69 kDa.

ターゲット情報

機能	Transcriptional regulator involved in development.
組織特異性	Expressed in thymus.
関連疾患	Defects in FOYN1 are the cause of T-cell immunodeficiency congenital alopecia and nail dystrophy (TIDAND) [MIM:601705]. A disorder characterized by the association of congenital alopecia, severe T-cell immunodeficiency, and ridging and pitting of all nails.
配列類似性	Contains 1 fork-head DNA-binding domain.
細胞内局在	Nucleus.

画像



Anti-FOYN1 antibody (ab96374) at 1/1000 dilution + Raji whole cell lysate at 30 µg

Predicted band size: 69 kDa

7.5% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors