

Anti-FSBA antibody [2F5/5] ab18575

製品の概要

製品名	Anti-FSBA antibody [2F5/5]
製品の詳細	Mouse monoclonal [2F5/5] to FSBA
由来種	Mouse
アプリケーション	適用あり: ELISA, IP
種交差性	交差種: Species independent
特記事項	<p>FSBA 2F5/5 has been shown to detect Protein Kinase C-β1 treated with FSBA.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
精製度	Immunogen affinity purified
一次抗体 備考	FSBA 2F5/5 has been shown to detect Protein Kinase C-β1 treated with FSBA.
ポリ/モノ	モノクローナル
クローン名	2F5/5
ミエローマ	Sp2/0-Ag14
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab18575の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
IP		Use at an assay dependent dilution.

ターゲット情報

関連性 FSBA, a reactive ATP analogue, can be used to selectively label protein kinases, which may subsequently be detected and affinity purified using an anti-FSBA antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors