

### Anti-GATA3 antibody ab106625

★★★★★ [1 Abreviews](#) [27 References](#) [画像数 1](#)

#### 製品の概要

製品名	Anti-GATA3 antibody
製品の詳細	Rabbit polyclonal to GATA3
由来種	Rabbit
特異性	At least two isoforms of GATA3 are known to exist; ab106625 will recognize both isoforms.
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide corresponding to Human GATA3 (C terminal). Database link: <a href="#">P23771</a>
ポジティブ・コントロール	WB: Human brain tissue lysate.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

## The Abpromise guarantee


**Abpromise保証は、次のテスト済みアプリケーションにおけるab106625の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 48 kDa.

## ターゲット情報

機能	Transcriptional activator which binds to the enhancer of the T-cell receptor alpha and delta genes. Binds to the consensus sequence 5'-AGATAG-3'.
組織特異性	T-cells and endothelial cells.
関連疾患	Defects in GATA3 are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia (HDR) [MIM:146255]; also known as Barakat syndrome.
配列類似性	Contains 2 GATA-type zinc fingers.
細胞内局在	Nucleus.
製品の状態	There are 2 isoforms produced by alternative splicing.

## 画像



**Lane 1 :** Anti-GATA3 antibody (ab106625) at 1 µg/ml

**Lane 2 :** Anti-GATA3 antibody (ab106625) at 2 µg/ml

**All lanes :** Human brain tissue lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size:** 48 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors