

### Anti-Gentamicin antibody ab123989

#### 製品の概要

---

製品名	Anti-Gentamicin antibody
製品の詳細	Sheep polyclonal to Gentamicin
由来種	Sheep
アプリケーション	<b>適用あり:</b> ELISA, Competitive ELISA
種交差性	<b>交差種:</b> Species independent
免疫原	Gentamycin conjugated to BSA.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 0.88% Sodium chloride, PBS
精製度	IgG fraction
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

---

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab123989の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

