

Anti-GIP antibody [04] ab30679

★★★★★ [4 Abreviews](#) [8 References](#) [画像数 2](#)

製品の概要

製品名	Anti-GIP antibody [04]
製品の詳細	Mouse monoclonal [04] to GIP
由来種	Mouse
特異性	Ab30679 recognises GIP. It cross reacts up to 10% with PACAP and GLP2 but <0.1% with glucagon, GLP1 or VIP.
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Synthetic human GIP coupled to carrier. The immunogen was adsorbed onto Al(OH) ₃ and emulsified in Freund's incomplete adjuvant.
特記事項	<p>Fusion partner: SP2/mIL6.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituents: 0.0268% PBS, 0.87% Sodium chloride</p>
精製度	Protein A purified
一次抗体 備考	Fusion partner: SP2/mIL6.
ポリ/モノ	モノクローナル
クローン名	04

アイソタイプ IgG1
軽鎖の種類 kappa

アプリケーション


The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab30679の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (1)	1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報


機能 Potent stimulator of insulin secretion and relatively poor inhibitor of gastric acid secretion.
配列類似性 Belongs to the glucagon family.
細胞内局在 Secreted.

画像



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-GIP antibody [04]
(ab30679)

ab30679 at 1/100 staining GIP in human duodenum by IHC-P.



ab30679 at 1/100 staining GIP in human small intestine by IHC-P.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-GIP antibody [04] (ab30679)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors