



Anti-GLUT5 antibody ab36057

9 References

製品の概要

製品名	Anti-GLUT5 antibody
製品の詳細	Rabbit polyclonal to GLUT5
由来種	Rabbit
特異性	Does not react with GLUT1, 2, 3 or 4 from human, mouse or rat.
アプリケーション	適用あり: IHC-Fr, IHC-P, WB
種交差性	交差種: Human 非交差種: Mouse, Rat
免疫原	Synthetic peptide corresponding to Human GLUT5 aa 490-500 (C terminal) conjugated to keyhole limpet haemocyanin. Sequence: ELKEPPVTSEQ
	 Run BLAST with  Run BLAST with
特記事項	Concentration is approximately 5 mg/ml. The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 0.05% BSA, PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab36057の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-Fr		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

ターゲット情報

機能	Cytochalasin B-sensitive carrier. Seems to function primarily as a fructose transporter.
組織特異性	Expressed in small intestine, and at much lower levels in kidney, skeletal muscle, and adipose tissue.
配列類似性	Belongs to the major facilitator superfamily. Sugar transporter (TC 2.A.1.1) family. Glucose transporter subfamily.
細胞内局在	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors