

Anti-Glycophorin C/GPC antibody [Ret40f] ab9521

画像数 2

製品の概要

製品名	Anti-Glycophorin C/GPC antibody [Ret40f]
製品の詳細	Mouse monoclonal [Ret40f] to Glycophorin C/GPC
由来種	Mouse
特異性	This antibody is specific to beta sialoglycophorin also known as Glycophorin C/GPC.
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Tissue, cells or virus corresponding to Glycophorin C/GPC. Red cell ghosts. Database link: P04921
特記事項	<p>Glycophorin C/GPC is found on the human erythrocyte membrane and is one of four sialic acid rich polypeptide.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.30 Preservative: 0.05% Sodium azide
精製度	Tissue culture supernatant
一次抗体 備考	Glycophorin C is found on the human erythrocyte membrane and is one of four sialic acid rich polypeptide.
ポリ/モノ	モノクローナル
クローン名	Ret40f

ミエローマ	unknown
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab9521の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/15 - 1/30. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

ターゲット情報


機能	This protein is a minor sialoglycoprotein in human erythrocyte membranes. The blood group Gerbich antigens and receptors for Plasmodium falciparum merozoites are most likely located within the extracellular domain. Glycophorin-C plays an important role in regulating the stability of red cells.
組織特異性	Glycophorin-C is expressed in erythrocytes. Glycophorin-D is ubiquitous.
配列類似性	Belongs to the glycophorin-C family.
細胞内局在	Cell membrane. Linked to the membrane via band 4.1.

画像



Formalin-fixed, paraffin-embedded human tonsil tissue stained for Glycophorin C/GPC using ab9521 at 1/50 dilution in immunohistochemical analysis. Antigen retrieval with citrate buffer pH 6.0.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Glycophorin C/GPC antibody [Ret40f] (ab9521)



Immunohistochemical analysis of formalin-fixed, paraffin-embedded Human tonsil tissue, staining Glycophorin C/GPC with ab9521.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Glycophorin C/GPC antibody [Ret40f] (ab9521)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors