

Anti-GNPDA2 antibody ab106363

[2 References](#) [画像数 2](#)

製品の概要

製品名	Anti-GNPDA2 antibody
製品の詳細	Rabbit polyclonal to GNPDA2
由来種	Rabbit
特異性	ab106363 is predicted not to cross-react with GNPDA1.
アプリケーション	適用あり: WB, ICC
種交差性	交差種: Human
免疫原	An 8 amino acid synthetic peptide near the C terminal of Human GNPDA2 (NP_612208).
ポジティブ・コントロール	WB: Human kidney lysate; ICC: Human GNPDA2 cells.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab106363の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 31 kDa.
ICC		Use a concentration of 4 µg/ml.

ターゲット情報

配列類似性 Belongs to the glucosamine/galactosamine-6-phosphate isomerase family.
細胞内局在 Cytoplasm.

画像



Anti-GNPDA2 antibody (ab106363) at 1 µg/ml + Human kidney lysate at 15 µg

Predicted band size: 31 kDa

Western blot - Anti-GNPDA2 antibody (ab106363)



Immunocytochemical analysis of Human GNPDA2 in 293 cells stained with ab106363 at 4 µg/mL.

Immunocytochemistry - Anti-GNPDA2 antibody (ab106363)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors