

Gold Anti-hCG beta antibody [M705159] ab31206

製品の概要

製品名	Gold Anti-hCG beta antibody [M705159]
製品の詳細	Gold Mouse monoclonal [M705159] to hCG beta
由来種	Mouse
標識	Gold
アプリケーション	適用あり: ELISA
種交差性	交差種: Human
免疫原	Full length native protein (purified). This information is proprietary to Abcam and/or its suppliers.
特記事項	Labeled with 40nm colloidal gold particles gold (~10ug/ml).

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.4 Preservative: 0.09% Sodium azide Constituents: PBS, 1% BSA
精製度	Affinity purified
特記事項(精製)	Purified by immunoaffinity chromatography prior to labeling.
ポリクローナル	モノクローナル
クローン名	M705159
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab31206の使用に適用されず**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.

ターゲット情報

関連性 Function: Stimulates the ovaries to synthesize the steroids that are essential for the maintenance of pregnancy. Tissue specificity: Placenta. Similarity: Belongs to the glycoprotein hormones subunit beta family. Developmental stage: Made by the first trimester placenta.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors