


Anti-GORASP2/GRASP55 antibody ab74579

2 References [画像数 1](#)

製品の概要

製品名	Anti-GORASP2/GRASP55 antibody
製品の詳細	Rabbit polyclonal to GORASP2/GRASP55
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat 
免疫原	Synthetic peptide corresponding to Human GORASP2/GRASP55 (internal sequence).
ポジティブ・コントロール	HeLa, HepG2, Jurkat and LOVO cell extracts.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: 50% Glycerol, 0.87% Sodium chloride, PBS</p>
精製度	Without Mg ²⁺ and Ca ²⁺
特記事項 (精製)	Immunogen affinity purified
ポリ/モノ	ab74579 was affinity purified from rabbit antiserum by affinity chromatography using epitope specific immunogen.
アイソタイプ	ポリクローナル
	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab74579の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Predicted molecular weight: 47 kDa.

ターゲット情報

機能	May be involved in assembly and membrane stacking of the Golgi cisternae, and in the process by which Golgi stacks reform after mitotic breakdown. May regulate the intracellular transport and presentation of a defined set of transmembrane proteins, such as transmembrane TGFA.
配列類似性	Belongs to the GORASP family. Contains 1 PDZ (DHR) domain.
翻訳後修飾	Myristoylated. Myristoylation is essential for the Golgi targeting. Palmitoylated. Phosphorylated in mitotic cells.
細胞内局在	Golgi apparatus membrane.

画像



All lanes : Anti-GORASP2/GRASP55 antibody (ab74579) at 1/500 dilution

Lane 1 : HeLa cell extracts

Lane 2 : HepG2 cell extracts

Lane 3 : Jurkat cell extracts

Lane 4 : LOVO cell extracts

Lane 5 : HeLa cell extracts with immunizing peptide at 10 µg

Lysates/proteins at 10 µg per lane.

Predicted band size: 47 kDa

Observed band size: 47 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors