


Anti-GPx-7 antibody ab96257

面格粒

医薬用外劇物

製品の概要

製品名	Anti-GPx-7 antibody
製品の詳細	Rabbit polyclonal to GPx-7
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Cow, Monkey 
免疫原	Recombinant fragment corresponding to Human GPx-7 aa 1-186. Database link: NP_056511
ポジティブ・コントロール	GPx-7 transfected 293T cell lysate; H441
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.21% Tris, 0.75% Glycine, 10% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab96257の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。


アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 21 kDa.
IHC-P		1/100.

ターゲット情報

配列類似性 Belongs to the glutathione peroxidase family.

細胞内局在 Secreted.

画像




All lanes : Anti-GPx-7 antibody (ab96257) at 1/500 dilution

Lane 1 : Non transfected 293T cell lysate

Lane 2 : GPx-7 transfected 293T cell lysate

Predicted band size: 21 kDa

12% SDS-PAGE



Immunohistochemical analysis of paraffin-embedded H441, using ab96257 antibody at 1/100 dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors