

Anti-Granzyme A antibody [GA6] ab10870

★★★★★ [1 Abreviews](#) [4 References](#)

製品の概要

製品名	Anti-Granzyme A antibody [GA6]
製品の詳細	Mouse monoclonal [GA6] to Granzyme A
由来種	Mouse
アプリケーション	適用あり: IHC-P, IHC-FoFr, WB
種交差性	交差種: Human
免疫原	Recombinant full length protein corresponding to Human Granzyme A. Database link: P12544
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Protein A purified
特記事項(精製)	Purified IgG prepared by affinity chromatography from tissue culture supernatant.
ポリ/モノ	モノクローナル
クローン名	GA6
ミエローマ	Sp2/0
アイソタイプ	IgG1

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab10870の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★ (1)	Use at an assay dependent concentration.
IHC-FoFr		Use at an assay dependent concentration. PubMed: 19827750
WB		Use at an assay dependent concentration. Predicted molecular weight: 29 kDa. Recognises a 60kD disulphide-linked homodimeric protein.

ターゲット情報

機能	This enzyme is necessary for target cell lysis in cell-mediated immune responses. It cleaves after Lys or Arg. Cleaves APEX1 after 'Lys-31' and destroys its oxidative repair activity. Involved in apoptosis.
配列類似性	Belongs to the peptidase S1 family. Granzyme subfamily. Contains 1 peptidase S1 domain.
細胞内局在	Secreted. Cytoplasmic granule.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors