



Anti-HA tag antibody [12CA5] ab1424

★★★★★ [12 Abreviews](#) [59 References](#) [画像数 4](#)

製品の概要

製品名	Anti-HA tag antibody [12CA5]
製品の詳細	Mouse monoclonal [12CA5] to HA tag
由来種	Mouse
特異性	The antibody detects HA-tag fused to the C-terminus or N-terminus of targeted proteins in transfected cells. The antibody works in ChIP but less well than ab9110 (see ChIP data).
アプリケーション	適用あり: ICC/IF, ChIP, IP, WB
種交差性	交差種: Species independent
免疫原	Synthetic peptide. Derived from the human influenza hemagglutinin protein Database link: ADG21028
	 Run BLAST with  Run BLAST with
ポジティブ・コントロール	WB: Rab6-HA transfected HEK-293T whole cell lysate. IP: Rab6-HA transfected HEK-293T whole cell lysate. ICC/IF: Baby hamster kidney cell expressing HA-tag CFTR protein. ChIP: Xenopus laevis oocytes injected with mRNA for HA-tagged human BORIS.
特記事項	<p>This product is manufactured by BioVision, an Abcam company and was previously called 3996 HA-Tag Antibody (Clone 12CA5[2]).</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Constituent: 1% Dibasic monohydrogen potassium phosphate
キャリア・フリー	はい
精製度	Affinity purified

ポリ/モノ	モノクローナル
クローン名	12CA5
アイソタイプ	IgG2b
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab1424の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
ICC/IF	★★★★★ (1)	Use a concentration of 10 µg/ml.
ChIP		Use a concentration of 3 µg/ml. Use 3 µg per 25 µg chromatin (see figure legend).
IP	★★★★★ (5)	Use a concentration of 10 µg/ml.
WB	★★★★★ (6)	Use a concentration of 1 µg/ml.

ターゲット情報

関連性

Human influenza hemagglutinin (HA) is a surface glycoprotein required for the infectivity of the human virus. The HA tag is derived from the HA molecule corresponding to amino acids 98-106 has been extensively used as a general epitope tag in expression vectors. Many recombinant proteins have been engineered to express the HA tag, which does not appear to interfere with the bioactivity or the biodistribution of the recombinant protein. This tag facilitates the detection, isolation, and purification of the proteins.

画像



Anti-HA tag antibody [12CA5] (ab1424) at 1/1000 dilution +
HEK293T whole cell lysate over-expressing HA-tagged Rab6 at 20
µg

Secondary


HRP-conjugated goat anti-mouse polyclonal IgG at 1/5000 dilution

Developed using the ECL technique.


Performed under reducing conditions.

Exposure time: 5 seconds

Blocked with 5% milk for 30 min at 25°C



This picture shows baby hamster kidney cell expressing HA-tag CFTR protein. This type of cells were grown on coverslip and immunofluorescence was performed. Dilution of 1:200 was used and incubated for 2hrs with antibody.




Immunoprecipitation - Anti-HA tag antibody [12CA5]
(ab1424)

Data courtesy of an anonymous AbReview

Recombinant HA-tagged Rab6 was immunoprecipitated from 1 mg transfected HEK-293T (human epithelial cell line from embryonic kidney transformed with large T antigen) whole cell extract with 2 μ L ab1424.

Lane 1: Rab6-HA transfected HEK-293T whole cell extract (Input)

Lane 2: ab1424 IP in Rab6-HA transfected HEK-293T whole cell extract.



ChIP - Anti-HA tag antibody [12CA5] (ab1424)

Xenopus laevis oocytes were injected with mRNA for HA-tagged human BORIS. Chromatin was prepared according to the Abcam X-ChIP protocol. Oocytes were fixed with formaldehyde for 10min. The ChIP was performed with 25 μ g of chromatin, 3 μ g of ab1424 (anti-HA, light blue) and 3 μ g of **ab18337** (anti-BORIS, dark blue), and 20 μ L of Protein A/G sepharose beads. A non-specific antibody was used as a control (yellow). The immunoprecipitated DNA was quantified by real time PCR (Taqman approach).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors