

Anti-Hairy antibody [1/24] ab20165

★★★★★ **1 Abreviews** **3 References**

製品の概要

| | |
|----------|---|
| 製品名 | Anti-Hairy antibody [1/24] |
| 製品の詳細 | Mouse monoclonal [1/24] to Hairy |
| 由来種 | Mouse |
| アプリケーション | 適用あり: ChIP, IHC, WB |
| 種交差性 | 交差種: Drosophila melanogaster |
| 免疫原 | Other Immunogen Type corresponding to Hairy. |
| 特記事項 | <p>This antibody stains Drosophila embryos.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|---------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | Preservative: 0.02% Sodium azide Constituent: 99.98% PBS |
| 精製度 | Protein A/G purified |
| 一次抗体 備考 | This antibody stains Drosophila embryos |
| ポリ/モノ | モノクローナル |
| クローン名 | 1/24 |
| ミエローマ | Sp2/0-Ag14 |
| アイソタイプ | IgG1 |
| 軽鎖の種類 | unknown |

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab20165の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| ChIP | | Use at an assay dependent concentration. |
| IHC | | Use at an assay dependent concentration. |
| WB | | Use at an assay dependent concentration. Predicted molecular weight: 37 kDa. |

ターゲット情報

| | |
|-------|--|
| 関連性 | Pair rule protein that regulates embryonic segmentation and adult bristle patterning. Transcriptional repressor of genes that require a bHLH protein for their transcription (e.g. the Fushi tarazu gene. Transcription repression requires formation of a complex with a corepressor |
| 細胞内局在 | Nuclear |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors