

Anti-HCV subtype 1b NS5B antibody ab35586

★★★★★ [1 Abreviews](#) [17 References](#) [画像数 1](#)

製品の概要

| | |
|----------|---|
| 製品名 | Anti-HCV subtype 1b NS5B antibody |
| 製品の詳細 | Rabbit polyclonal to HCV subtype 1b NS5B |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: ELISA, WB |
| 種交差性 | 交差種: Hepatitis C virus |
| 免疫原 | Recombinant fragment (Hepatitis C Virus) with a 21 amino acid C-terminal truncation |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | Constituent: PBS |
| 精製度 | Protein A purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab35586の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| ELISA | ★★★★★ (1) | Use at an assay dependent dilution. |
| WB | | 1/2000. Detects a band of approximately 64 kDa. |

ターゲット情報

細胞内局在

Host cytoplasm, nucleus, mitochondrion, endoplasmatic reticulum; secreted.

画像



All lanes : Anti-HCV subtype 1b NS5B antibody (ab35586) at 1/2000 dilution

Lane 1 : Lysate from Huh7 cells, mock transfection

Lane 2 : Lysate from Huh7 cells transfected with NS5B

Observed band size: 64 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors