

# Anti-HGF antibody [EGH2 4C12.1] ab10829

**1 References** [画像数 1](#)

### 製品の概要

製品名	Anti-HGF antibody [EGH2 4C12.1]
製品の詳細	Mouse monoclonal [EGH2 4C12.1] to HGF
由来種	Mouse
アプリケーション	<b>適用あり:</b> Sandwich ELISA
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment corresponding to Human HGF.
特記事項	<p>Plasminogen related growth factors are potent effectors of growth, movement and differentiation of epithelia and endothelia. Membrane receptor for HGF/SF is the product of the proto-oncogene c-met a tyrosine kinase.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
精製度	Protein A purified
一次抗体 備考	Plasminogen related growth factors are potent effectors of growth, movement and differentiation of epithelia and endothelia. Membrane receptor for HGF/SF is the product of the proto-oncogene c-met a tyrosine kinase.
ポリ/モノ	モノクローナル

クローン名	EGH2 4C12.1
ミエローマ	NS0
アイソタイプ	IgG2b

## アプリケーション


**The Abpromise guarantee** **Abpromise保証は、**次のテスト済みアプリケーションにおけるab10829の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
<b>Sandwich ELISA</b>		Use a concentration of 5 µg/ml. Can be paired for Sandwich ELISA with <b>Rabbit polyclonal to HGF (ab24865)</b> . For sandwich ELISA, use this antibody as Capture at 5µg/ml with <b>ab24865</b> as Detection.

## ターゲット情報

<b>機能</b>	Potent mitogen for mature parenchymal hepatocyte cells, seems to be a hepatotrophic factor, and acts as a growth factor for a broad spectrum of tissues and cell types. Activating ligand for the receptor tyrosine kinase MET by binding to it and promoting its dimerization.
<b>関連疾患</b>	Deafness, autosomal recessive, 39
<b>配列類似性</b>	Belongs to the peptidase S1 family. Plasminogen subfamily. Contains 4 kringle domains. Contains 1 PAN domain. Contains 1 peptidase S1 domain.

## 画像



Standard Curve for HGF (Analyte: **ab632**) dilution range 1pg/ml to 1ug/ml using Capture Antibody Mouse monoclonal [EGH2 4C12.1] to HGF (ab10829) at 5ug/ml and Detector Antibody Rabbit polyclonal to HGF (**ab24865**) at 0.5ug/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors