

### Anti-HIPK2 antibody [EPR3819] ab108543

リコンビナント **RabMAb**

**6 References**   [画像数 4](#)

#### 製品の概要

製品名	Anti-HIPK2 antibody [EPR3819]
製品の詳細	Rabbit monoclonal [EPR3819] to HIPK2
由来種	Rabbit
特異性	The immunogen used for this product shares 69% homology and 6 consecutive amino acids with HIPK1. Cross-reactivity with this protein has not been confirmed experimentally.
アプリケーション	<b>適用あり:</b> Flow Cyt (Intra), WB <b>適用なし:</b> IHC-P
種交差性	<b>交差種:</b> Mouse, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	HT-1080 and HepG2 whole cell lysate ( <b>ab7900</b> )
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> <li>- High batch-to-batch consistency and reproducibility</li> <li>- Improved sensitivity and specificity</li> <li>- Long-term security of supply</li> <li>- Animal-free production</li> </ul> <p>For more information <a href="#">see here</a>.</p> <p>Our RabMAb<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb<sup>®</sup> patents</a>.</p> <p>Rat: We have preliminary internal testing data to indicate this antibody may not react with this species. Please contact us for more information.</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituents: 0.1% BSA, 40% Glycerol (glycerin, glycerine), 9.85% Tris glycine, 50% Tissue culture supernatant</p>

精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR3819
アイソタイプ	IgG

## アプリケーション

**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab108543の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。


アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		1/1000 - 1/10000. <b>ab172730</b> - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
WB		1/2000 - 1/10000. Detects a band of approximately 100, 130 kDa (predicted molecular weight: 131 kDa).

**追加情報** Is unsuitable for IHC-P.

## ターゲット情報

機能	Protein kinase acting as a corepressor of several transcription factors, including SMAD1 and POU4F1/Bm3a and probably NK homeodomain transcription factors. Inhibits cell growth and promotes apoptosis. Involved in transcriptional activation of TP53 and TP73. Phosphorylation of TP53 may be mediated by a TP53-HIPK2-AXIN1 complex. In response to TGFB, cooperates with DAXX to activate JNK. Phosphorylates the antiapoptotic factor CTBP1 and promotes its proteasomal degradation. In the Wnt/beta-catenin signaling pathway acts as an intermediate kinase between TAK1 and NLK to promote the proteasomal degradation of MYB (By similarity). Phosphorylates CBX4 upon DNA damage and promotes its E3 SUMO-protein ligase activity. PML, HIPK2 and FBXO3 may act synergically to activate p53/TP53-dependent transactivation.
組織特異性	Highly expressed in heart, muscle and kidney. Weakly expressed in a ubiquitous way. Down-regulated in several thyroid and breast tumors.
配列類似性	Belongs to the protein kinase superfamily. CMGC Ser/Thr protein kinase family. HIPK subfamily. Contains 1 protein kinase domain.
翻訳後修飾	Phosphorylated on tyrosines (By similarity). Autophosphorylated. Sumoylated. When conjugated it is directed to nuclear speckles. Desumoylated by SENP1 (By similarity). Sumoylation on Lys-32 is promoted by the E3 SUMO-protein ligase CBX4. Ubiquitinated by FBXO3, leading to rapid proteasome-dependent degradation. This degradation, but not ubiquitination, is prevented in the presence of PML.
細胞内局在	Nucleus > PML body. Cytoplasm. Concentrated in PML/POD/ND10 nuclear bodies. Small amounts are cytoplasmic.

## 画像



Western blot - Anti-HIPK2 antibody [EPR3819]  
(ab108543)

Anti-HIPK2 antibody [EPR3819] (ab108543) at 1/2000 dilution +  
Mouse heart tissue lysate at 20 µg

### Secondary


Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/10000 dilution

**Predicted band size:** 131 kDa

**Observed band size:** 100,130 kDa


**Exposure time:** 180 seconds

**Blocking and diluting buffer and concentration:** 5% NFDM  
/TBST.



Flow Cytometry (Intracellular) - Anti-HIPK2 antibody  
[EPR3819] (ab108543)

Overlay histogram showing HEK293 cells stained with ab108543 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab108543, 1/10000 dilution) for 30 min at 22°C. The secondary antibody used was Alexa Fluor® 488 goat anti-rabbit IgG (H&L) ([ab150077](#)) at 1/2000 dilution for 30 min at 22°C. Isotype control antibody (black line) was rabbit IgG (monoclonal) (0.1µg/1x10<sup>6</sup> cells) used under the same conditions. Unlabelled sample (blue line) was also used as a control. Acquisition of >5,000 events were collected using a 20mW Argon ion laser (488nm) and 525/30 bandpass filter. This antibody gave a positive signal in HEK293 cells fixed with 4% paraformaldehyde (10 min)/permeabilized with 0.1% PBS-Tween for 20 min used under the same conditions.



Western blot - Anti-HIPK2 antibody [EPR3819]  
(ab108543)

**All lanes :** Anti-HIPK2 antibody [EPR3819] (ab108543) at 1/2000 dilution

**Lane 1 :** HT-1080 cell lysate

**Lane 2 :** HepG2 cell lysate

Lysates/proteins at 10 µg per lane.

### Secondary





**All lanes :** HRP labelled goat anti-rabbit at 1/2000 dilution

**Predicted band size:** 131 kDa

**Observed band size:** 130 kDa

**Additional bands at:** 100 kDa (possible isoform)

Why choose a recombinant antibody?

 <p><b>Research with confidence</b> Consistent and reproducible results</p>	 <p><b>Long-term and scalable supply</b> Recombinant technology</p>
 <p><b>Success from the first experiment</b> Confirmed specificity</p>	 <p><b>Ethical standards compliant</b> Animal-free production</p>

Anti-HIPK2 antibody [EPR3819] (ab108543)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors