




Anti-Histone H3 (acetyl K9) antibody - ChIP Grade ab4441

★★★★★ [46 Abreviews](#) [415 References](#) [画像数 4](#)

製品の概要

製品名	Anti-Histone H3 (acetyl K9) antibody - ChIP Grade
製品の詳細	Rabbit polyclonal to Histone H3 (acetyl K9) - ChIP Grade
由来種	Rabbit
特異性	In Dot blot detects 50ng of mono-acetylated peptide corresponding to position Lys9 in the N-terminal sequence of Histone H3. Does not detect the mono-acetylated peptide corresponding to acetyl-lysine at position 14 or unacetylated Histone H3.
アプリケーション	適用あり: WB, ChIP
種交差性	交差種: Human 交差が予測される動物種: <i>Saccharomyces cerevisiae</i> 
免疫原	Synthetic peptide corresponding to Human Histone H3 aa 1-12 (acetyl K9) conjugated to Keyhole Limpet Haemocyanin (KLH). Sequence: ARTKQTAR (Ac) KSTG-C <div>  Run BLAST with  Run BLAST with </div>
ポジティブ・コントロール	ChIP: U2OS cells, Human Cell lysate - nuclear (Primary CD34+ cells) WB: Mouse MEL cell nuclear lysate, HeLa cell acid-extract
特記事項	<p>Learn about ChIP assay kits, other ChIP antibodies, protocols and more in the ChIP assay guide.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

バッファー	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 0.184% Tris glycine, 30% Glycerol, 0.87% Sodium chloride
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab4441の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (20)	1/10000. Detects a band of approximately 17 kDa.
ChIP	★★★★★ (12)	Use 2µg for 10 ⁶ cells. Use GAPDH primer pair ab267832 as positive control.

ターゲット情報


機能	Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.
配列類似性	Belongs to the histone H3 family.
発生段階	Expressed during S phase, then expression strongly decreases as cell division slows down during the process of differentiation.
翻訳後修飾	<p>Acetylation is generally linked to gene activation. Acetylation on Lys-10 (H3K9ac) impairs methylation at Arg-9 (H3R8me2s). Acetylation on Lys-19 (H3K18ac) and Lys-24 (H3K24ac) favors methylation at Arg-18 (H3R17me).</p> <p>Citrullination at Arg-9 (H3R8ci) and/or Arg-18 (H3R17ci) by PAD4 impairs methylation and represses transcription.</p> <p>Asymmetric dimethylation at Arg-18 (H3R17me2a) by CARM1 is linked to gene activation. Symmetric dimethylation at Arg-9 (H3R8me2s) by PRMT5 is linked to gene repression. Asymmetric dimethylation at Arg-3 (H3R2me2a) by PRMT6 is linked to gene repression and is mutually exclusive with H3 Lys-5 methylation (H3K4me2 and H3K4me3). H3R2me2a is present at the 3' of genes regardless of their transcription state and is enriched on inactive promoters, while it is absent on active promoters.</p> <p>Methylation at Lys-5 (H3K4me), Lys-37 (H3K36me) and Lys-80 (H3K79me) are linked to gene activation. Methylation at Lys-5 (H3K4me) facilitates subsequent acetylation of H3 and H4. Methylation at Lys-80 (H3K79me) is associated with DNA double-strand break (DSB) responses and is a specific target for TP53BP1. Methylation at Lys-10 (H3K9me) and Lys-28 (H3K27me) are linked to gene repression. Methylation at Lys-10 (H3K9me) is a specific target for HP1 proteins (CBX1, CBX3 and CBX5) and prevents subsequent phosphorylation at Ser-11 (H3S10ph) and acetylation of H3 and H4. Methylation at Lys-5 (H3K4me) and Lys-80 (H3K79me)</p>

require preliminary monoubiquitination of H2B at 'Lys-120'. Methylation at Lys-10 (H3K9me) and Lys-28 (H3K27me) are enriched in inactive X chromosome chromatin. Phosphorylated at Thr-4 (H3T3ph) by GSG2/haspin during prophase and dephosphorylated during anaphase. Phosphorylation at Ser-11 (H3S10ph) by AURKB is crucial for chromosome condensation and cell-cycle progression during mitosis and meiosis. In addition phosphorylation at Ser-11 (H3S10ph) by RPS6KA4 and RPS6KA5 is important during interphase because it enables the transcription of genes following external stimulation, like mitogens, stress, growth factors or UV irradiation and result in the activation of genes, such as c-fos and c-jun. Phosphorylation at Ser-11 (H3S10ph), which is linked to gene activation, prevents methylation at Lys-10 (H3K9me) but facilitates acetylation of H3 and H4. Phosphorylation at Ser-11 (H3S10ph) by AURKB mediates the dissociation of HP1 proteins (CBX1, CBX3 and CBX5) from heterochromatin. Phosphorylation at Ser-11 (H3S10ph) is also an essential regulatory mechanism for neoplastic cell transformation. Phosphorylated at Ser-29 (H3S28ph) by MLTK isoform 1, RPS6KA5 or AURKB during mitosis or upon ultraviolet B irradiation. Phosphorylation at Thr-7 (H3T6ph) by PRKCBB is a specific tag for epigenetic transcriptional activation that prevents demethylation of Lys-5 (H3K4me) by LSD1/KDM1A. At centromeres, specifically phosphorylated at Thr-12 (H3T11ph) from prophase to early anaphase, by DAPK3 and PKN1. Phosphorylation at Thr-12 (H3T11ph) by PKN1 is a specific tag for epigenetic transcriptional activation that promotes demethylation of Lys-10 (H3K9me) by KDM4C/JMJD2C. Phosphorylation at Tyr-42 (H3Y41ph) by JAK2 promotes exclusion of CBX5 (HP1 alpha) from chromatin. Monoubiquitinated by RAG1 in lymphoid cells, monoubiquitination is required for V(D)J recombination (By similarity). Ubiquitinated by the CUL4-DDB-RBX1 complex in response to ultraviolet irradiation. This may weaken the interaction between histones and DNA and facilitate DNA accessibility to repair proteins.


細胞内局在

Nucleus. Chromosome.

画像



Chromatin was prepared from U2OS cells according to the Abcam X-ChIP protocol. Cells were fixed with formaldehyde for 10 min. The ChIP was performed with 25 µg of chromatin, 2 µg of **ab4441** (blue), and 20 µl of protein A/G sepharose beads. No antibody was added to the beads control (yellow). The immunoprecipitated DNA was quantified by real time PCR (Taqman approach). Primers and probes are located in the first kb of the transcribed region.



ChIP analysis using unpurified ab4441 binding Histone H3 in human primary CD34+ cell lysate. Cells were cross-linked for 15 minutes with 1% formaldehyde. Samples were incubated with primary antibody (0.2 µg/µg of chromatin) for 16 hours at 4°C. Protein binding was detected using real-time PCR.


Positive controls:

Region upstream of the transcription start site of the ACTB gene.

Region upstream of the transcription start site of the EGR1 gene.

Negative Controls:

Region on chromosome 12 (Untr12) that is far from any known gene annotation and not expected to be bound by Histone H3 (acetyl K9).



All lanes : Anti-Histone H3 (acetyl K9) antibody - ChIP Grade (ab4441) at 1/2500 dilution

Lane 1 : Mouse MEL cell nuclear lysate

Lane 2 : Mouse MEL cell nuclear lysate treated with 0.4 µM Trichostatin A treatment for 18 hours.

Lysates/proteins at 9 µg per lane.


Secondary

All lanes : Donkey Anti-Rabbit IgG H&L (HRP) ([ab6802](#)) at 1/20000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Exposure time: 1 minute



All lanes : Anti-Histone H3 (acetyl K9) antibody - ChIP Grade (ab4441) at 0.5 µg/ml

Lane 1 : Untreated HeLa cell acid-extract

Lane 2 : HeLa cell acid-extract treated with sodium butyrate

ab4441, Histone H3, acetylated (Lys9) Pab Rabbit x Human

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors