


Anti-Histone H4 antibody [mAbcam 31830] - ChIP Grade ab31830

★★★★★ [6 Abreviews](#) [51 References](#) [画像数 7](#)

製品の概要

製品名	Anti-Histone H4 antibody [mAbcam 31830] - ChIP Grade
製品の詳細	Mouse monoclonal [mAbcam 31830] to Histone H4 - ChIP Grade
由来種	Mouse
アプリケーション	適用あり: Flow Cyt (Intra), IHC-P, WB, ChIP, ICC/IF 適用なし: IP
種交差性	交差種: Cow, Human 交差が予測される動物種: Mouse, Rat 
免疫原	Synthetic peptide corresponding to Human Histone H4 aa 50 to the C-terminus (C terminal) conjugated to keyhole limpet haemocyanin. Database link: P62805 (Peptide available as ab13843)
ポジティブ・コントロール	This antibody gave a positive signal in the following lysates: Calf Thymus Histone Preparation Nuclear Lysate HeLa Histone Preparation Nuclear Lysate Histone H4 Recombinant Protein IHC-P: FFPE human breast fibroadenoma. ICC/IF: HeLa cell line
特記事項	<p>This antibody clone is manufactured by Abcam. If you require a custom buffer formulation or conjugation for your experiments, please contact orders@abcam.com.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

バッファー	pH: 7.50 Preservative: 0.02% Sodium azide Constituents: PBS, 6.97% L-Arginine
精製度	IgG fraction
ポリ/モノ	モノクローナル
クローン名	mAbcam 31830
ミエローマ	Sp2/0-Ag14
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab31830の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		Use 1 µg for 10 ⁶ cells. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.
IHC-P	★★★★★ (1)	Use a concentration of 0.05 - 1 µg/ml. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.
WB	★★★★★ (4)	Use a concentration of 1 µg/ml. Detects a band of approximately 13 kDa (predicted molecular weight: 14 kDa). Can be blocked with Human Histone H4 peptide (ab13843) .
ChIP		Use 5 µg for 25 µg of chromatin.
ICC/IF		Use a concentration of 5 µg/ml.

追加情報 Is unsuitable for IP.

ターゲット情報

機能	Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.
配列類似性	Belongs to the histone H4 family.
翻訳後修飾	Acetylation at Lys-6 (H4K5ac), Lys-9 (H4K8ac), Lys-13 (H4K12ac) and Lys-17 (H4K16ac) occurs in coding regions of the genome but not in heterochromatin. Citruination at Arg-4 (H4R3ci) by PADI4 impairs methylation. Monomethylation and asymmetric dimethylation at Arg-4 (H4R3me1 and H4R3me2a, respectively) by PRMT1 favors acetylation at Lys-9 (H4K8ac) and Lys-13 (H4K12ac). Demethylation is performed by JMJD6. Symmetric dimethylation on Arg-4 (H4R3me2s) by the

PRDM1/PRMT5 complex may play a crucial role in the germ-cell lineage.


Monomethylated, dimethylated or trimethylated at Lys-21 (H4K20me1, H4K20me2, H4K20me3). Monomethylation is performed by SET8. Trimethylation is performed by SUV420H1 and SUV420H2 and induces gene silencing.

Ubiquitinated by the CUL4-DDB-RBX1 complex in response to ultraviolet irradiation. This may weaken the interaction between histones and DNA and facilitate DNA accessibility to repair proteins. Monoubiquitinated at Lys-92 of histone H4 (H4K91ub1) in response to DNA damage. The exact role of H4K91ub1 in DNA damage response is still unclear but it may function as a licensing signal for additional histone H4 post-translational modifications such as H4 Lys-21 methylation (H4K20me).

Sumoylated, which is associated with transcriptional repression.

細胞内局在 Nucleus. Chromosome.

画像




Immunocytochemistry/ Immunofluorescence - Anti-Histone H4 antibody [mAbcam 31830] - ChIP Grade (ab31830)

ab31830 staining Histone H4 in HeLa cells. The cells were fixed with 100% methanol (5 min), permeabilized with 0.1% PBS-Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1%PBS-Tween for 1h. The cells were then incubated overnight at 4°C with ab31830 at 1µg/ml and **ab6046**, Rabbit polyclonal to beta Tubulin - Loading Control. Cells were then incubated with **ab150117**, Goat polyclonal Secondary Antibody to Mouse IgG H&L (Alexa Fluor® 488) preadsorbed at 1/1000 dilution (shown in green) and **ab150080**, Goat polyclonal Secondary Antibody to Rabbit IgG - H&L (Alexa Fluor® 594) at 1/1000 dilution (shown in pseudocolour magenta). Nuclear DNA was labelled with DAPI (shown in blue).

Also suitable in cells fixed with 4% paraformaldehyde (10 min).

Image was acquired with a high-content analyser (Operetta CLS, Perkin Elmer) and a maximum intensity projection of confocal sections is shown.



All lanes : Anti-Histone H4 antibody [mAbcam 31830] - ChIP Grade (ab31830) at 1 µg/ml

Lane 1 : Calf Thymus Histone Preparation Nuclear Lysate at 0.5 µg/ml

Lane 2 : HeLa Nuclear Lysate (Triton enriched) at 10 µg/ml

Lane 3 : Histone H4 Recombinant Protein at 0.1 µg/ml

Lane 4 : Histone H3.1 Recombinant Protein at 0.1 µg/ml

Secondary

All lanes : Goat polyclonal to Mouse IgG - H&L - Pre-Adsorbed (HRP) at 1/5000 dilution


Developed using the ECL technique.

Performed under reducing conditions.


Predicted band size: 14 kDa

Additional bands at: 13 kDa. We are unsure as to the identity of these extra bands.

Exposure time: 10 seconds



Chromatin was prepared from HeLa cells according to the Abcam X-ChIP protocol. Cells were fixed with formaldehyde for 10 minutes. The ChIP was performed with 25µg of chromatin, 5µg of ab31830 (blue), and 20µl of protein A/G sepharose beads. No antibody was added to the beads control (yellow). The immunoprecipitated DNA was quantified by real time PCR (Taqman approach for active and inactive loci, Sybr green approach for heterochromatic loci). Primers and probes are located in the first kb of the transcribed region.



Flow Cytometry (Intracellular) - Anti-Histone H4 antibody [mAbcam 31830] - ChIP Grade (ab31830)

Overlay histogram showing HeLa cells stained with ab31830 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab31830, 1µg/1x10⁶ cells) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) ([ab96879](#)) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG1 [ICIGG1] ([ab91353](#), 2µg/1x10⁶ cells) used under the same conditions. Acquisition of >5,000 events was performed.



Western blot - Anti-Histone H4 antibody [mAbcam 31830] - ChIP Grade (ab31830)

All lanes : Anti-Histone H4 antibody [mAbcam 31830] - ChIP Grade (ab31830) at 1 µg/ml

Lane 1 : Calf Thymus Histone Preparation Nuclear Lysate at 0.5 µg

Lane 2 : HeLa Histone Preparation Nuclear Lysate at 2.5 µg

Lane 3 : Histone H4 Recombinant Protein at 0.1 µg

Lane 4 : Histone H3.1 Recombinant Protein at 0.1 µg


Secondary

All lanes : Goat polyclonal to Mouse IgG - H&L - Pre-Adsorbed (HRP) at 1/3000 dilution

Performed under reducing conditions.

Predicted band size: 14 kDa

Observed band size: 13 kDa



All lanes : Anti-Histone H4 antibody [mAbcam 31830] - ChIP Grade (ab31830) at 5 µg/ml

Lane 1 : Calf Thymus Histone Preparation Nuclear Lysate at 0.5 µg

Lane 2 : HeLa Histone Preparation Nuclear Lysate at 2.5 µg

Lane 3 : Histone H4 Recombinant Protein at 0.1 µg

Lane 4 : Histone H3.1 Recombinant Protein at 0.1 µg

Lane 5 : Histone H2A Recombinant Protein at 0.1 µg

Lane 6 : Histone H2B Recombinant Protein at 0.1 µg

Lane 7 : Histone H1 Recombinant Protein at 0.1 µg

Lane 8 : Calf Thymus Histone Preparation Nuclear Lysate at 0.5 µg with Human Histone H4 peptide ([ab13843](#)) at 5 µg/ml

Lane 9 : HeLa Histone Preparation Nuclear Lysate at 2.5 µg with Human Histone H4 peptide ([ab13843](#)) at 5 µg/ml

Lane 10 : Histone H4 Recombinant Protein at 0.1 µg with Human Histone H4 peptide ([ab13843](#)) at 5 µg/ml

Lane 11 : Histone H3.1 Recombinant Protein at 0.1 µg with Human Histone H4 peptide ([ab13843](#)) at 5 µg/ml

Lane 12 : Histone H2A Recombinant Protein at 0.1 µg with Human Histone H4 peptide ([ab13843](#)) at 5 µg/ml

Lane 13 : Histone H2B Recombinant Protein at 0.1 µg with Human Histone H4 peptide ([ab13843](#)) at 5 µg/ml

Lane 14 : Histone H1 Recombinant Protein at 0.1 µg with Human Histone H4 peptide ([ab13843](#)) at 5 µg/ml

Secondary

All lanes : Goat polyclonal to Mouse IgG - H&L - Pre-Adsorbed (HRP) at 1/3000 dilution

Performed under reducing conditions.

Predicted band size: 14 kDa

Observed band size: 13 kDa

Exposure time: 2 minutes



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Histone H4 antibody [mAbcam 31830] - ChIP Grade (ab31830)

IHC image of Histone H4 staining in human breast fibroadenoma formalin fixed paraffin embedded tissue section, performed on a Leica Bond™ system using the standard protocol F. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH6, epitope retrieval solution 1) for 20 mins. The section was then incubated with [ab31830](#), 0.1µg/ml, for 15 mins at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors