

Anti-HLA-DQ antibody [B-K27] - BSA and Azide free ab47342

画像数 1

製品の概要

製品名	Anti-HLA-DQ antibody [B-K27] - BSA and Azide free
製品の詳細	Mouse monoclonal [B-K27] to HLA-DQ - BSA and Azide free
由来種	Mouse
特異性	ab47342 recognises both a 28 kDa and 33 kDa complex.
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	Tissue, cells or virus corresponding to Human HLA-DQ. T4+ CTL clone HG-38
ポジティブ・コントロール	Human lymphocytes
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.30 Constituent: 100% PBS
キャリア・フリー	Sterile-filtered through 0.22 µm. Carrier and preservative free はい
精製度	Ion Exchange Chromatography
ポリ/モノ	モノクローナル
クローン名	B-K27
ミエローマ	x63-Ag8.653

アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション


The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab47342の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration. ab170191 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

ターゲット情報

関連性	Human Leukocyte Antigens are highly polymorphic proteins that are involved in the presentation of antigens to the T-cell receptor. There are two classes of HLA antigens: class I (HLA A, HLA B and HLA C) and class II (HLA D).
細胞内局在	Type I membrane protein

画像



Flow Cytometry analysis of lymphocytes labeling HLA-DQ with ab47342.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors