

Anti-HOXA7 antibody [743C1a] ab51235

2 References [画像数 1](#)

製品の概要

| | |
|----------|---|
| 製品名 | Anti-HOXA7 antibody [743C1a] |
| 製品の詳細 | Mouse monoclonal [743C1a] to HOXA7 |
| 由来種 | Mouse |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Recombinant fragment |
| 免疫原 | Recombinant fragment corresponding to Human HOXA7. |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|-----------|---|
| 製品の状態 | Liquid |
| 保存方法 | Store at +4°C short term (1-2 weeks) (add antibacterial agent). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | <p>pH: 7.40</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituents: 1% BSA, PBS</p> |
| 精製度 | Protein G purified |
| 特記事項 (精製) | ab51235 was purified from culture supernatant of the hybridoma, cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum. |
| ポリ/モノ | モノクローナル |
| クローン名 | 743C1a |
| アイソタイプ | IgG1 |

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab51235の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| WB | | Use at an assay dependent concentration. Predicted molecular weight: 25 kDa. |

ターゲット情報

| | |
|-------|---|
| 機能 | Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. |
| 配列類似性 | Belongs to the Antp homeobox family. Contains 1 homeobox DNA-binding domain. |
| 細胞内局在 | Nucleus. |

画像



Anti-HOXA7 antibody [743C1a] (ab51235)

Predicted band size: 25 kDa

Observed band size: 31 kDa

Western blot analysis of immunised recombinant protein

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors