

### Anti-HOXB4 antibody [EP1919Y] ab76093

リコンビナント **RabMAb**

**2 References**    **画像数 8**

#### 製品の概要

製品名	Anti-HOXB4 antibody [EP1919Y]
製品の詳細	Rabbit monoclonal [EP1919Y] to HOXB4
由来種	Rabbit
アプリケーション	<b>適用あり:</b> Flow Cyt (Intra), WB, IP, IHC-P, ICC/IF
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide within Human HOXB4 aa 100-200 (internal sequence). The exact sequence is proprietary.
ポジティブ・コントロール	WB: K-562 lysates. IHC-P: human lung carcinoma tissue. ICC/IF: K-562 cells. Flow Cyt (intra): K-562 cells. IP: K-562 cells.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> <li>- High batch-to-batch consistency and reproducibility</li> <li>- Improved sensitivity and specificity</li> <li>- Long-term security of supply</li> <li>- Animal-free production</li> </ul> <p>For more information <a href="#">see here</a>.</p> <p>Our RabMAb<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb<sup>®</sup> patents</a>.</p> <p>Mouse, Rat: We have preliminary internal testing data to indicate this antibody may not react with these species. Please contact us for more information.</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA</p>
精製度	Protein A purified

ポリ/モノ	モノクローナル
クローン名	EP1919Y
アイソタイプ	IgG

## アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab76093の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		1/100. <b>ab172730</b> - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
WB		1/1000 - 1/10000. Predicted molecular weight: 28 kDa.
IP		1/60. <b>For unpurified use please use at 1/40</b>
IHC-P		1/100 - 1/250. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. See <b><u>IHC antigen retrieval protocols</u></b>
ICC/IF		1/100 - 1/250.

## ターゲット情報


機能	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
配列類似性	Belongs to the Antp homeobox family. Deformed subfamily. Contains 1 homeobox DNA-binding domain.
発生段階	Expressed in whole embryos and fetuses at 5-9 weeks from conception.
細胞内局在	Nucleus.

## 画像




Immunocytochemistry/ Immunofluorescence - Anti-HOXB4 antibody [EP1919Y] (ab76093)

Immunocytochemistry/ Immunofluorescence analysis of K-562 (Human chronic myelogenous leukemia lymphoblast) cells labeling HOXB4 with purified ab76093 at 1:100 dilution (10 µg/ml). Cells were fixed in 4% Paraformaldehyde and permeabilized with 0.1% tritonX-100. Cells were counterstained with Ab195889 Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) 1:200 (2.5 µg/ml). Goat anti rabbit IgG (Alexa Fluor® 488, **ab150077**) was used as the secondary antibody at 1:1000 (2 µg/ml) dilution. DAPI (blue) was used as nuclear counterstain. PBS instead of the primary antibody was used as the secondary antibody only control.




Immunoprecipitation - Anti-HOXB4 antibody [EP1919Y] (ab76093)

ab76093 (Purified) at 1:60 dilution (2 µg) immunoprecipitating HOXB4 in K-562 whole cell lysate.  
 Lane 1 (input): K-562 (Human chronic myelogenous leukemia lymphoblast) whole cell lysate 10 µg  
 Lane 2 (+): ab76093 & K-562 whole cell lysate  
 Lane 3 (-): Rabbit monoclonal IgG (**ab172730**) instead of ab76093 in K-562 whole cell lysate  
 For western blotting, VeriBlot for IP Detection Reagent (HRP) (**ab131366**) was used at 1:1000 dilution.  
 Blocking and diluting buffer: 5% NFDm/TBST.




Flow Cytometry (Intracellular) - Anti-HOXB4 antibody  
[EP1919Y] (ab76093)

Intracellular Flow Cytometry analysis of K-562 (Human chronic myelogenous leukemia lymphoblast) cells labeling HOXB4 with purified ab76093 at 1/100 dilution (10 µg/ml) (Red). Cells were fixed with 4% Paraformaldehyde and permeabilised with 90% Methanol. A Goat anti rabbit IgG (Alexa Fluor® 488, [ab150077](#)) secondary antibody was used at 1/2000. Isotype control - Rabbit monoclonal IgG (Black). Unlabeled control - Cell without incubation with primary antibody and secondary antibody (Blue).



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-HOXB4 antibody  
[EP1919Y] (ab76093)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human lung carcinoma tissue sections labeling HOXB4 with purified ab76093 at 1:100 dilution (11.86 µg/ml). Heat mediated antigen retrieval was performed using [ab93684](#) (Tris/EDTA buffer, pH 9.0). ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain.



Western blot - Anti-HOXB4 antibody [EP1919Y]  
(ab76093)


Anti-HOXB4 antibody [EP1919Y] (ab76093) at 1/1000 dilution  
(Purified) + K-562 (Human chronic myelogenous leukemia  
lymphoblast) whole cell lysates at 15 µg

#### Secondary

Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution


**Predicted band size:** 28 kDa

**Observed band size:** 34 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-  
embedded sections) - Anti-HOXB4 antibody  
[EP1919Y] (ab76093)





Immunohistochemical analysis of paraffin-embedded human  
placental tissue, using unpurified ab76093, at 1/100 dilution.  
  
Perform heat mediated antigen retrieval before commencing with  
IHC staining protocol.



Flow Cytometry (Intracellular) - Anti-HOXB4 antibody  
[EP1919Y] (ab76093)

Overlay histogram showing K562 cells stained with unpurified ab76093 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab76093, 1/100 dilution) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-rabbit IgG (H+L) (**ab96899**) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was rabbit IgG (monoclonal) (1µg/1x10<sup>6</sup> cells) used under the same conditions. Acquisition of >5,000 events was performed.

Why choose a recombinant antibody?

 <p><b>Research with confidence</b> Consistent and reproducible results</p>	 <p><b>Long-term and scalable supply</b> Recombinant technology</p>
 <p><b>Success from the first experiment</b> Confirmed specificity</p>	 <p><b>Ethical standards compliant</b> Animal-free production</p>

Anti-HOXB4 antibody [EP1919Y] (ab76093)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- 
- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors