

HRP Anti-Digoxigenin antibody [HRP.21H8] ab6212

★★★★★ 4 Abbreviations 22 References 画像数 1

医薬用外劇物

製品の概要

製品名	HRP Anti-Digoxigenin antibody [HRP.21H8]
製品の詳細	HRP Mouse monoclonal [HRP.21H8] to Digoxigenin
由来種	Mouse
標識	HRP
アプリケーション	適用あり: IHC-Fr
種交差性	交差種: Species independent
免疫原	Chemical/ Small Molecule conjugated to bovine serum albumin.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.01% Thimerosal (merthiolate)</p> <p>Constituents: PBS, 50% Glycerol, 1% BSA</p>
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	HRP.21H8
ミエローマ	unknown
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab6212の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。


アプリケーション	Abreviews	特記事項
IHC-Fr	★★★★★ (2)	1/200 - 1/2000. 1/200 - 1/2000. Use in combination with tyramide signal amplification.

ターゲット情報

関連性

Digoxigenin (DIG) is a steroid found exclusively in the flowers and leaves of the plants *Digitalis purpurea* and *Digitalis lanata*. Digoxigenin is chemically closely related to Digoxin, the cardiac glycoside used for the treatment of various heart diseases. The term 'genin' at the end of Digoxigenin, refers to only the aglycone portion (without the sugar) part of the molecule, thus Digoxigenin is the steroid component of Digoxin, - minus the (digitose) sugar residues. DIG can be covalently added to proteins or nucleic acids which makes it very useful in diverse applications.

画像



ab6212 at a dilution of 1/1000 staining Digoxigenin from human pancreatic tissue sections by immunohistochemistry (frozen). The antibody was incubated with the tissue for 2 hours.

This image is courtesy of an Abreview submitted by an anonymous researcher on **9 December 2005**.

Immunohistochemistry (Frozen sections) - HRP
Anti-Digoxigenin antibody [HRP.21H8] (ab6212)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors