

HRP Anti-HIV1 gp120 antibody ab53840

★★★★★ 1 Abreviews 1 References

医薬用外劇物

製品の概要

製品名	HRP Anti-HIV1 gp120 antibody
製品の詳細	HRP Goat polyclonal to HIV1 gp120
由来種	Goat
標識	HRP
特異性	Gp120 is the predominant band in WB.
アプリケーション	適用あり: WB, ELISA, IHC-Fr
種交差性	交差種: Human immunodeficiency virus
免疫原	Full length native protein, from strain IIIB.
特記事項	<p>Estimated molar HRP:IgG substitution is 2-3. Free enzyme is absent.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.002% Thimerosal (merthiolate) Constituents: 1% BSA, PBS
特記事項(精製)	>95%
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab53840の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	1/20 - 1/200.
ELISA		1/200 - 1/1000.
IHC-Fr		1/20 - 1/200.

ターゲット情報

関連性

HIV1 is equipped with the envelope gp160 glycoprotein for interaction with Langerhans cells (LCs) and dendritic cells (DCs), the members of the innate immune system, which confront the virus at the portal of virus entry in the human body. These cells are equipped with receptors by which they bind and endocytose the virus. The gp120 glycoprotein is used for binding to CD4 receptor and CCR5 co-receptor of T helper 2 (Th2) cells, and is able to induce FcepsilonRI(+) hematopoietic cells to produce IL4, which inactivates the host adaptive immune response.

細胞内局在

Envelope surface glycoprotein.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors