

HRP Anti-Human IgA antibody [1H9] ab7383

7 References

製品の概要

製品名	HRP Anti-Human IgA antibody [1H9]
製品の詳細	HRP Mouse monoclonal [1H9] to Human IgA
由来種	Mouse
標識	HRP
特異性	IgA Fc region.
アプリケーション	適用あり: ELISA, Sandwich ELISA
種交差性	交差種: Human
免疫原	IgA from human colostrum.
特記事項	<p>Concentration varies from lot to lot and can be provided on request.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.05% Proclin 300</p> <p>Constituents: 0.8% Sodium chloride, 0.02% Potassium chloride, 0.18% Dibasic monohydrogen sodium phosphate, 0.024% Monobasic dihydrogen potassium phosphate</p>
精製度	Protein G purified
特記事項 (精製)	Purified by chromatography on protein G Sepharose. Purity is tested by electrophoresis.
ポリ/モノ	モノクローナル
クローン名	1H9
ミエローマ	Sp2/0

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab7383の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
Sandwich ELISA		Use at an assay dependent concentration. Can be paired for Sandwich ELISA with <u>Mouse monoclonal [3B7] to Human IgA (ab7400)</u> . Can be used as Detection antibody with recommended pair(s) as capture <u>ab7400</u> .

ターゲット情報

関連性	Ig alpha is the major immunoglobulin class in body secretions. It may serve both to defend against local infection and to prevent access of foreign antigens to the general immunologic system.
細胞内局在	Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors