

### HRP Anti-pCrebtide antibody ab21625

#### 製品の概要

---

製品名	HRP Anti-pCrebtide antibody
製品の詳細	HRP Rabbit polyclonal to pCrebtide
由来種	Rabbit
標識	HRP
特異性	Recognize pCrebtide only. Does not cross-react with unphosphorylated Crebtide
アプリケーション	<b>適用あり:</b> ELISA
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide corresponding to pCrebtide conjugated to keyhole limpet haemocyanin.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6 Preservative: 0.1% Sodium azide Constituent: Tris buffered saline
精製度	Immunogen affinity purified
特記事項(精製)	Affinity purified with the pCrebtide on Agarose and immunoaffinity adsorption with Crebtide on Agarose.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

---

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab21625の使用に適用されず**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.

## ターゲット情報

---

**関連性**      Crebtide is a peptide derived from the protein Creb, and it is well documented substrate for PKA and PKC kinases. pCrebtide is the phosphorylated product of these kinases.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors