

### Anti-Hsp40 antibody ab69402

★★★★★ [5 Abreviews](#) [13 References](#) [画像数 1](#)

#### 製品の概要

製品名	Anti-Hsp40 antibody
製品の詳細	Rabbit polyclonal to Hsp40
由来種	Rabbit
特異性	This antibody detects a ~40 kDa protein, corresponding to the apparent molecular mass of Hsp40 on SDS-PAGE immunoblots. It does not detect DnaJ in E.coli lysates.
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Mouse, Rat, Human, Recombinant fragment <b>非交差種:</b> Escherichia coli
免疫原	Recombinant human Hsp40
ポジティブ・コントロール	WB: Heat shocked HeLa cells and 3T3 cells, PC-12 cells.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 50% Glycerol (glycerin, glycerine), PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション


**The Abpromise guarantee**      Abpromise保証は、次のテスト済みアプリケーションにおけるab69402の使用に適用されます  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (2)	1/10000. Detects a band of approximately 40 kDa (predicted molecular weight: 38 kDa).

## ターゲット情報

機能	Interacts with HSP70 and can stimulate its ATPase activity. Stimulates the association between HSC70 and HIP.
配列類似性	Contains 1 J domain.
細胞内局在	Cytoplasm. Nucleus. Nucleus > nucleolus. Translocates rapidly from the cytoplasm to the nucleus, and especially to the nucleoli, upon heat shock.

## 画像



**All lanes :** Anti-Hsp40 antibody (ab69402) at 1/1000 dilution

**Lane 1 :** Molecular weight marker

**Lane 2 :** Hsp40 recombinant protein

**Lane 3 :** DnaJ recombinant protein

**Lane 4 :** Cell lysates prepared from heat shocked hela cells

**Lane 5 :** Cell lysates prepared from heat shocked 3T3 cells

**Lane 6 :** Cell lysates prepared from PC-12 cells

**Predicted band size:** 38 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors