

Anti-HSV1 gE Envelope Protein antibody [9H3] ab6510

9 References

製品の概要

製品名	Anti-HSV1 gE Envelope Protein antibody [9H3]
製品の詳細	Mouse monoclonal [9H3] to HSV1 gE Envelope Protein
由来種	Mouse
アプリケーション	適用あり: ELISA, ICC/IF, WB
種交差性	交差種: Herpes simplex virus
免疫原	Tissue, cells or virus corresponding to HSV1 gE Envelope Protein.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.4
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	9H3
ミエローマ	NS1/1-Ag4-1
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab6510の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration. 0.100 @ >1:102,400
ICC/IF		1/6400. 1+ @ 1:6,400
WB		Use at an assay dependent concentration. PubMed: 16571817

ターゲット情報

関連性 There are 2 antigenic types of Herpes Simplex Virus, HSV1 and HSV2 which share antigenic cross reactivity but different neutralization patterns and tend to produce different clinical symptoms. Humans are believed to be the natural host for HSV, but the virus is also capable of infecting various animals, including rodents.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors