

# Anti-HSV2 gG Envelope antibody [H1206] ab53470

## 1 References

### 製品の概要

---

製品名	Anti-HSV2 gG Envelope antibody [H1206]
製品の詳細	Mouse monoclonal [H1206] to HSV2 gG Envelope
由来種	Mouse
特異性	ab53470 is reactive with gG of Herpes Simplex Virus 2.
アプリケーション	<b>適用あり:</b> ICC/IF, WB
種交差性	<b>交差種:</b> Human herpesvirus 2
免疫原	Full length protein corresponding to HSV2 gG Envelope .
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.40 Constituent: PBS
精製度	Protein G purified
特記事項(精製)	Purified from TCS
ポリモノ	モノクローナル
クローン名	H1206
アイソタイプ	IgG1
軽鎖の種類	kappa

## アプリケーション

---

**The Abpromise guarantee** Abpromise保証は、 次のテスト済みアプリケーションにおけるab53470の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 10 µg/ml.
WB		Use a concentration of 10 µg/ml.

## ターゲット情報

---

**関連性** HSV2 is a member of the herpes virus family, Herpesviridae. All herpes viruses are morphologically identical. They have a large double stranded DNA genome. The virion consists of an icosahedral nucleocapsid which is surrounded by a lipid bilayer envelope. Between the capsid and the envelope is an amorphous layer of proteins, termed the tegument. Glycoprotein G (gG) is one of seven envelope glycoproteins in HSV2.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors