

# Anti-HURP antibody ab84509

**3 References** [画像数 3](#)

### 製品の概要

製品名	Anti-HURP antibody
製品の詳細	Rabbit polyclonal to HURP
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IHC-P, WB
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide corresponding to Human HURP aa 200-300. Database link: <a href="#">Q15398</a>
ポジティブ・コントロール	IHC-P: Human basal cell carcinoma and human seminoma tissues.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.1% BSA, Tris buffered saline
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

The Abpromise guarantee

**Abpromise保証は、** 次のテスト済みアプリケーションにおけるab84509の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/100 - 1/500.
WB		1/200.

## ターゲット情報


機能	Potential cell cycle regulator that may play a role in carcinogenesis of cancer cells. Mitotic phosphoprotein regulated by the ubiquitin-proteasome pathway. Key regulator of adherens junction integrity and differentiation that may be involved in CDH1-mediated adhesion and signaling in epithelial cells.
組織特異性	Abundantly expressed in fetal liver. Expressed at lower levels in bone marrow, testis, colon, and placenta.
配列類似性	Belongs to the SAPAP family.
発生段階	Elevated levels of expression detected in the G2/M phase of synchronized cultures of HeLa cells.
翻訳後修飾	Ubiquitinated, leading to its degradation. Phosphorylated upon DNA damage, probably by ATM or ATR. Decreased phosphorylation levels are associated with the differentiation of intestinal epithelial cells.
細胞内局在	Nucleus. Cytoplasm. Cytoplasm > cytoskeleton > spindle. Localizes to the spindle in mitotic cells. Colocalizes with CDH1 at sites of cell-cell contact in intestinal epithelial cells.

## 画像




ab84509, at 1/250 dilution, staining HURP in formalin-fixed, paraffin-embedded human basal cell carcinoma by Immunohistochemistry, using DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-HURP antibody (ab84509)



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-HURP antibody (ab84509)

ab84509, at 1/250 dilution, staining HURP in formalin-fixed, paraffin-embedded human seminoma by Immunohistochemistry, using DAB staining.



Western blot - Anti-HURP antibody (ab84509)

Anti-HURP antibody (ab84509) at 1/200 dilution + 25µg of WCE from U2OS cells synchronized by thymidine arrest and release to enrich in mitotic cells (same on lane 1 and 2).

*Image kindly provided by Patrizia Lavia, CNR National Research Council, University "La Sapienza", Rome, Italy.*

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors